XENON
ready

TANGENS 2C

ANÄSTHESIESTATION
ANAESTHESIA UNIT

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WE DEVELOP WHAT YOU HAVE WISHED FOR

We have listened well. The experience of anaesthesia teams, their requirements and requests regarding a safe and reliable appliance therefore affect the design of the TANGENS 2C. And the operating efficiency of this Xenon-device oblige the challenges of cost pressure in the health care system. For the safety of your patients and efficiency in procedures.

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FLEXIBLE COMPACT SOLUTION

Integrating a patient monitor in a standard device is rather unusual so far. The compact TANGENS 2C offers you this advantage. With the aid of a folding arm this one can be put in any position for operating and surveying tasks. Moreover the compact construction of this anaesthetic unit considerably increases your space of movement.

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THE FUTURE OF ANAESTHESIA WITH XENON HAS STARTED

For the first time economic solution available

As precondition for a cost-conscious application of Xenon as anaesthetic EKU uses the electronic EGAMIX-system in its TANGENS 2C. The precise and extremely tight electronic gas mixer guarantees patient-compatible dosing and avoids expensive gas losses. Furthermore the integrated neurological monitoring system ANATREND for supervising the depth of anaesthesia contributes to optimal dosing of Xenon. Due to its reliable supervision the TANGENS 2C controls a precise supply also of other anaesthetics and thus ensures an individual anaesthesia management.

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A REAL EYE-CATCHER: THE CENTRAL SCREEN

Highest functionality and easy handling

The TANGENS 2C features integrated patient monitoring. With complete alarm control of the vital parameters. The operation is dialogue driven via six operating keys only. For eventually you neither have the time to read many-sided operating manuals nor to enter all standards. The TANGENS 2C executes standard routines for you in the background and grants more time for your patients due to its self-explanatory operation system.

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FLEXIBLE SOLUTIONS FOR YOUR ANAESTHESIA WORKSTATIONS

Innovative anaesthesia information system

EKU offers you a complete hard- and software solution for continuous documentation of the anaesthesia records. Including network connection to the clinical PDMS-system and support in the service diagnosis. The corresponding control panels can be mounted directly at the TANGENS-trolley or installed as a stand-alone-solution.

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USER-FRIENDLY UP TO THE SMALLEST DETAIL

The circle system of the TANGENS 2C is easy to use and to clean. The exchange of the absorber canister is done comfortably and fast. The system remains sealed even during the exchange.

The connections for the integrated patient monitoring are clearly arranged, functional and well accessible in the trolley's profile.

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HIGHER QUALITY IN ANAESTHESIA

Xenon as inhalation anaesthetics considerably expands the spectrum of the currently common methods for maintaining a general anaesthesia. The combination of inhaling and intravenous anaesthetics is supervised all together by the TANGENS 2C in its central screen.

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CEILING VERSION OF TANGENS 2C DVE

Due to its narrow design and its light self-weight the TANGENS 2C DVE is suited for installation in medium load ceiling-mounted supply systems. Directly connected to the ceiling-mounted supply system even the previously extensive docking apparatus do not apply anymore - another advantage for an ergonomic configuration of smaller operating theatres.

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OPTIMAL OPERATION INCREASES PRODUCTIVITY

We have enquired what counts in practice. This way we receive useful solutions such as the lockable drawer unit and the keyboard lighting by illuminated LEDs at the bottom frame of the central screen.

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General Data:	
dimensions:	600 x 730 x 1650 mm (W x D x H)
weight:	approx. 100 kg
display:	coloured, TFT-LCD, 12" or 15"
Power Supply:	
voltage:	100–240V ~ 50/60 Hz
power consumption:	max. 80VA (depending on the configuration)
internal accumulator:	12V/12Ah, approx. 2 hours battery operation
Circle System IBS:	
dimensions:	320 x 210 x 535 mm (W x D x H)
weight:	base: 3,0 kg IBS-body: 4,0 kg absorber: 1,2 kg
i-SORB CO ₂ absorber:	volume: 1420 ml, capacity: approx. 900 g
APL-Valve:	adjustments: SP, 5 hPa to 75 hPa
accuracy:	± 7 hPa at 4 L/min
Ventilator:	
Parameters:	
tidal volume (CMV):	20–1500 mL (flow: 3–160 L/min)
inspiratory flow (PCV):	10–80 L/min
frequency:	6–60 1/min
SIMV-frequency:	3–60 1/min
I:E-ratio:	1:0,5–1:4 (Ti: 0,3–10 s)
inspiratory pause:	0–50% Ti
PEEP:	OFF, 3–20 hPa
PMAX:	(PEEP +5)–80 hPa
flow trigger (spontaneous mode):	0,3–15 L/min
ventilation modes:	CMV, PCV, SIMV MMV, ASB in preparation
Accuracy:	
tidal volume:	± 10% (output „bellows in bottle“)
PMAX:	± 2 hPa or 5% rel.
PEEP:	± 1 hPa
driving gas:	
type:	compressed air or oxygen
pressure:	2–10 bar
flow:	max. 160 L/min
consumption:	approx. equal to the respective minute volume
Gas Mixer (normal mode):	
Parameters:	
gases:	O ₂ , AIR, N ₂ O
pressure:	3,5–6,0 bar
flow:	9,9 L/min
Fresh Gas:	
blend:	O ₂ /AIR or O ₂ /N ₂ O
flow (standby):	0 L/min
flow (operation):	0,2–9,9 L/min
resolution flow:	0,1 L/min
accuracy flow:	± 10% rel.
O ₂ -concentration:	25–100%
resolution O ₂ :	1 Vol%
accuracy O ₂ :	± 10% rel.
Settings:	- total fresh gas flow - O ₂ -concentration - blend (O ₂ /AIR or O ₂ /N ₂ O)
Security System:	
	- O ₂ -parameter limits - N ₂ O-blockade - O ₂ -emergency dosing - failure gas supply (failure AIR: 100% O ₂ , failure N ₂ O: 100% O ₂ , failure O ₂ : 100% AIR)
Gas Mixer (Xenon mode):	
Parameters:	
gases:	O ₂ , AIR, Xe
pressure:	3,5–6,0 bar
flow:	9,9 L/min
Fresh Gas:	
blend:	O ₂ /AIR or O ₂ /Xe
flow (standby):	0 L/min
flow (operation):	0,2–9,9 L/min
resolution flow:	0,01 L/min
accuracy flow:	± 10% rel.
O ₂ -concentration:	25–100% (30–100% in automatic mode)
resolution O ₂ :	1 Vol%
accuracy O ₂ :	± 10% rel.
Settings:	- blend (O ₂ /AIR or O ₂ /Xe) - normal mode: - total fresh gas flow - O ₂ -concentration

- automatic mode:	- total fresh gas flow - O ₂ -concentration
Security System:	
	- O ₂ -parameter limits - Xe-blockade - O ₂ -emergency dosing - failure gas supply (failure AIR: 100% O ₂ , failure Xe: O ₂ /AIR, failure O ₂ : 100% AIR)
Monitoring (normal mode):	
Anaesthetic Agents:	
anaesthetic agents:	Halothane, Enflurane, Isoflurane, Sevoflurane, Desflurane (optional: automatic agent identification)
measuring ranges:	HAL, ISO: 0–8,5 Vol%, ENF, SEV: 0–10 Vol%, DES: 0–22 Vol%
accuracy:	± (0,15 Vol% + 15% rel.)
resolution:	0,1 Vol%
MAC value display:	minimal alveolar concentration
N₂O:	
measuring range:	0–100 Vol%
accuracy:	± (2 Vol% + 8% rel.)
resolution:	1 Vol%
CO₂:	
measuring range:	0–10 Vol%
accuracy:	± 0,5 Vol% or ± 12% rel.
resolution:	0,1 Vol% (display also in mmHg)
Oxygen measurement paramagnetic:	
only available in connection with a module for automatic anaesthetic agent identification	
measuring range:	0–100 Vol%
accuracy:	± 1 Vol% or ± 2% rel.
resolution:	1 Vol%
Oxygen measurement by Galvanic Cell:	
measuring range:	0–100 Vol%
accuracy:	± 1 Vol% or ± 2% rel.
resolution:	1 Vol%
T90:	approx. 14 seconds
sensor life time:	approx. 1–4 years, depending on the measured oxygen concentration
Volume:	
Measured Values:	minute volume, tidal volume, breathing rate
Measuring Ranges:	
tidal volume:	0–2000 mL
minute volume:	0–99 L/min
breathing rate:	0–100 BPM
Accuracy:	
tidal volume:	± 8% or ± 0,1 L
minute volume:	± 8% or ± 0,5 L/min
breathing rate:	± 1 BPM
Pressure:	
measuring range:	-10–+90 hPa (display also in cmH ₂ O or mbar)
accuracy:	± 4% or ± 1hPa
resolution:	1 hPa
Monitoring (Xenon mode):	
Xe:	
measuring range:	0–100 Vol%
accuracy:	± 3 Vol%
resolution:	1 Vol%
CO₂:	
measuring range:	0–10 Vol%
accuracy:	± 0,5 Vol% or ± 12% rel.
resolution:	0,1 Vol% (display also in mmHg)
Oxygen Measurement:	
measuring range:	10–100 Vol%
accuracy:	± (1 Vol% + 2% rel.)
resolution:	1 Vol%
T90:	approx. 14 seconds
sensor life time:	approx. 1–4 years, depending on the measured oxygen concentration
Volume:	
Measured Values:	minute volume, tidal volume, breathing rate
Measuring Ranges:	
tidal volume:	0–2000 mL
minute volume:	0–99 L/min
breathing rate:	0–100 BPM
Accuracy:	
tidal volume:	± 15% or ± 0,1 L
minute volume:	± 15% or ± 0,5 L/min
breathing rate:	± 1 BPM

Above listed data are subject of change without further notice.

EKU Elektronik GmbH
Am Sportplatz
56291 Leiningen
Tel.: +49 (0) 6746-10-18
Fax: +49 (0) 6746-84-84
E-Mail: info@eku-elektronik.de
www.eku-elektronik.de